





Planmeca ClinicPatrol™






Water filtration system

Betriebsanleitung
Installation Instructions
Mode d'emploi
Istruzioni per l'uso
Gebruiksaanwijzing
Brugsanvisning
Bruksanvisningar
Käyttöohjeet

Copyright © 2024 BWT Holding GmbH.
Alle Rechte vorbehalten. All rights reserved. Tous droits réservés. Tutti i diritti riservati. Alle rechten voorbehouden. Alle rettigheder forbeholdes.

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Table of contents

1	General information	19
1.1	Abbreviations and subject index	19
1.2	Scope of delivery.....	19
1.3	Explanation of warnings.....	20
1.4	Symbols on the type plate.....	20
2	Technical data	20
2.1	Dimensions, connections and operating conditions for Planmeca ClinicPatrol	20
2.2	Operating conditions for Planmeca Patrol 14 membrane	21
2.3	Operating conditions for Planmeca Patrol 5	21
3	Use and function	21
3.1	Intended use.....	21
3.2	Design and function of the Planmeca ClinicPatrol	22
4	Operating and safety instructions	22
4.1	Responsibility of the operator	23
4.2	Warranty and exclusion of liability	23
4.3	Qualified personnel.....	23
4.4	Pressure	24
5	Installation and assembly	24
5.1	Requirements for installation	24
5.2	Unpacking the Planmeca ClinicPatrol	24
5.3	Hydraulic installation.....	24
6	Operating the reverse osmosis unit	25
6.1	Switching on the RO device.....	25
6.2	Hygiene concept and stoppages	25
6.3	Removing/replacing the RO cartridge	25
6.3.1	Disassembling/assembling the Planmeca Patrol 5	25
6.3.2	Disassembling/assembling the Planmeca Patrol 14 membrane	25
6.4	Installation and operation of the RO app.....	26
6.4.1	Installing the app	26
6.4.2	Using the app	26
7	Troubleshooting	26
7.1	Overview of the status and alert LED.....	26
7.2	Troubleshooting.....	27
8	Repairs and maintenance	29
8.1	Wearing parts	29
8.2	Cleaning	30
8.3	Standard IEC 60335-1	30
8.4	Disposal	30
9	Order numbers	30

1 General information

1.1 Abbreviations and subject index

Softener:	The water purification process removes the hardness from raw water. Hardness constituents are the portion of calcium and magnesium ions in the water.
Raw water:	The raw water is the untreated drinking water that is connected to the water inlet.
RO:	Abbreviation for reverse osmosis.
Permeate:	This is "pure water" that has been largely desalinated by reverse osmosis. The characteristic value is the electric conductivity in $\mu\text{S}/\text{cm}$.
Concentrate:	This is waste water containing the salts and minerals that have been removed from the raw water.
Membranes:	The "filter" of the device which is capable of desalinating the raw water by high pressure and flow.
TDS:	Total Dissolved Solids: Total content of dissolved salt, measured in mg/l .
SDI:	Silt Density Index: The "Silt Density Index" is a measure for the level of inlet water impurity.
Conductivity, electrical conductivity:	The smaller the value of the electrical conductivity measured (in $\mu\text{S}/\text{cm}$) by the RO device, the lower the salt concentration in the permeate product.
IOM:	Abbreviation for "Installation and Operating Manual"
Permeate yield (WCF):	The ratio between the amount of pure water (permeate) produced and the waste water produced. WCF stands for Water Conversion Factor.
Bypass setting with cold drinking water:	The bypass settings on the filter head have no function for this filter cartridge. The bypass setting of the bestdemin Plus must always be set to the 0 position for the device to function properly.

1.2 Scope of delivery

The following components are included in the scope of delivery of the RO device:

- Planmeca Patrol 14 membrane RO filter cartridge
- Planmeca Patrol 5 demineralisation
- Connection material: Connection cable for the power supply (type F / type I / type G), 3/4" x 10 mm adapter for plug-in connection, 3/8" x 8 mm adapter for plug-in connection, 8 mm concentrate facing for plug-in connection, 8 mm elbow connection piece for plug-in connection

The following components of the RO device can be seen on the front:

- 1 Planmeca Patrol 14 membrane RO filter cartridge
- 2 Planmeca Patrol 5 demineralisation
- 3 LED status
- 4 Top service cover

The following connections can be seen on the rear of the RO device:

- 5 Device ON/OFF switch
- 6 PE mains plug type IEC 320
- 7 Permeate connection 3/4"
- 8 Concentrate connection, 8 mm JG (5/16")
- 9 External tank connection 8 mm JG (5/16")
- 10 Permeate connection 3/8"
- 11 Technical data for the RO device

The top service cover can be removed to access the two integrated filter heads:

12 Filter head for the Planmeca Patrol 14 membrane with bypass setting, which must remain at position 0

13 Filter head for the Planmeca Patrol 5 without bypass setting

3

1.3 Explanation of warnings

⚠ DANGER!

Indicates a potentially dangerous situation that can lead to health impairments if not avoided.

⚠ WARNING!

Indicates a potentially dangerous situation that can lead to health impairments if not avoided.

⚠ CAUTION!

Indicates a potentially dangerous situation that can lead to property damage if not avoided.

① NOTE!

Highlights recommendations and information for efficient, trouble-free operation.

1.4 Symbols on the type plate

Pressure

Typical capacity

Temperature

Typical flow rate

Date of installation and replacement of the filter cartridge

2 Technical data

2.1 Dimensions, connections and operating conditions for Planmeca ClinicPatrol

Electrical connection / fuse protection	220–240 V / 50–60 Hz / 10 A
Internal device fuse	T1.25AL250V
Mains voltage fluctuations	Max. ± 10% of the rated voltage
Overvoltage category	II
Electrical power consumption (operation/standby)	200 W / < 3 W
Plug standard (grounded PE mains plug)	IEC-320 socket
Protection class	IP21
Water connection inlet	3/4" external thread
Water connection output (demineralised water)	3/8" external thread
Water connection for concentrate	8 mm connector
Water connection for external tank	8 mm connector
Dimensions: Width, depth, height (W×D×H)	277 mm x 297 mm x 505 mm
Weight, dry	19.0 kg
Weight, wet	22.6 kg
Usage	Suitable for indoor use only
Max. operating height	2000 m
Relative humidity	Max. 80% for T < 31 °C, decreasing linearly to 50% at T = 40 °C
Level of contamination	2

① NOTE!

The device is intended for intermittent operation with the following duty cycle: $T_{on} = 5 \text{ min} / T_{off} = 10 \text{ min}$ under typical conditions; the values for T_{on} and T_{off} may vary depending on the external conditions (ambient temperature, inlet water temperature, outlet pressure).

2.2 Operating conditions for Planmeca Patrol 14 membrane

Permeate production (amount produced) ^I	l/min (l/h)	2 (120) at 15 °C
Salt retention rate	%	> 97
Water conversion factor WCF (factory setting), ^{II} , ^{III}	%	50
Feedwater flow (input)	l/min (l/h)	min. 4 (240)
Nominal flow	l/h	120
Concentrate (waste outfeed)	l/min (l/h)	Approx. 2.0 (120)
Working pressure	bar	7
Intake water pressure	MPa (bar)	0.15 – 0.4 (1.5 – 4.0)
Water temperature (min./max.)	°C	4 – 30
Ambient temperature (min./max.)	°C	4 – 40
Iron + manganese (Fe+Mn)	mg/l	< 0.05
Silicate (SiO ₂)	mg/l	< 15
Salt content, total dissolved solids (TDS)	mg/l	< 500
Silt Density Index (SDI)	%/min	< 3
Oxidants	mg/l	< 0.05

^I The actual nominal flow rate may deviate slightly from the flow rate indicated in the table due to fluctuations in the feedwater quality, the flow pressure, the water temperature and the permeate counter-pressure (e.g. with greater permeate pumping heights).

^{II} As a general rule, the manufacturer recommends pre-treating the feedwater.

^{III} The RO device is factory-configured for a WCF of about 50%.

2.3 Operating conditions for Planmeca Patrol 5

Nominal flow	l/h	100
Working pressure range	bar	2 – 8
Intake water pressure	bar	> 1.2
Pressure loss at 30 l/h	bar	0.4
Pressure loss at 60 l/h	bar	0.6
Water temperature, min.–max.	°C	4 – 30
Ambient temperature, min.–max.	°C	4 – 40
Ambient temperature during storage/transport, min.–max.	°C	-20 – 40
Bed volumes	l	2.4
Operating position		Vertical
Conductivity limit of the permeate alarm	µS/cm	1 – 200
Demineralsised water conductivity at the outlet	µS/cm	< 3
Alarm at	µS/cm	1 – 100
Abort at	µS/cm	1 – 100

3 Use and function

3.1 Intended use

This RO device is a combination of reverse osmosis and a filter cartridge for Desalination of the permeate. The device

- may only be used for Desalination of cold water that meets the legal requirements for drinking water quality.
- produces desalinated water, which can be used as input water for autoclaves, steam sterilisation in laboratories and medicine and for electrolyzers for hydrogen production.
- filters particles and dissolved salts from the water.
- Filters particles out of the water.

Any other use is considered improper.

ⓘ NOTE!

The device is not suitable for continuous operation. There is a risk that the pump will overheat if the operating time is too long.

3.2 Design and function of the Planmeca ClinicPatrol

The RO device generates permeate which is desalinated in a second step with a mixed-bed ion exchanger. A simplified flow chart of the RO device is shown on the cover.

4

The water passes through the following stages:

- 1 Pump: Water is pumped into the RO filter cartridge at constant pressure.
- 2 Desalination: The water flows through the RO membrane
- 3 Permeate return: Permeate is recirculated until the conductivity is below an adjustable set value.
- 4 Ion exchanger: Demineralisation is carried out by the mixed-bed ion exchanger
- 5 System separator in accordance with DIN 1717 category 5
- 6 Transfer pump

4 Operating and safety instructions

Despite all safety precautions, residual risks remain with every product, especially if it is handled improperly. Every piece of technical equipment requires regular maintenance and repairs in order to function properly.

A basic precondition for safe working is the adherence to all stated safety and operating instructions. In addition, the local accident prevention provisions and the general safety provisions effective at the place of operation are applicable. The illustrations in these instructions are intended to provide a basic understanding and may deviate from the actual design of the device. Claims cannot be made based on this.

⚠ WARNING!

The RO device may not be operated with feedwater that is microbially contaminated or whose origin or quality is unknown.

Any deviation from proper use, such as desalination of feedwater of impermissible quality (non-drinking water), may result in irreversible damage to health or property (via microbial contamination of the RO device, for example).

Disconnect the RO device from the water supply prior to maintenance work on the drinking water supply. Flush the water line before reconnecting the RO device.

Before installation, disconnect the device and any terminal devices from the power (pull the mains plug).

If the device is not raised or carried correctly, there is a risk of injuries or damage to the device. Follow the instructions below:

- ▶ Consider the weight of the device before lifting
- ▶ Follow the regulations on protective clothing, such as safety shoes or anti-slip gloves
- ▶ If possible, reach under the device to carry it safely.
- ▶ If necessary, use a lifting or transport device.
- ▶ If necessary, ask other people to assist you.
- ▶ Secure the device during transport.

⚠ CAUTION!

Observe all national directives for drinking water installations (e.g. DIN 1988, EN 1717), general sanitary requirements and technical data for the protection of drinking water.

Improper installation of the RO device may cause damage to it.

The inflow water must not exceed the limit values given in the technical data or the calcium solubility limit!

The water fed into the device must be cold water that conforms to the statutory drinking water requirements and the quality requirements in *Technical data* [2].

Install a stop valve upstream of the RO device.

Only hoses that comply with DVGW W 543 may be used to connect the device.

If the product has been stored below 0°C, leave the unpacked product in the ambient temperature of the installation site for at least 24 hours before using it.

Do not install the filter system near heat sources and open flames.

The RO device may not come into contact with chemicals, solvents, or vapours.

The installation site must be free of frost and protected from direct sunlight.

① NOTE!

The materials were selected in accordance with the requirements of DIN 18879-1 and EN 14898.

The pressure resistance of the RO device corresponds to DIN 18879-1.

The desalinated water (permeate) produced must not be used for drinking.

Country-specific guidelines must be observed for the use of desalinated water (permeate).

Avoid leaving the device in storage for a long time in order to reduce the risk of idle contamination.

If the feedwater is treated with oxidising disinfectants (such as chlorine, chlorine dioxide, or similar), then an activated carbon filter must be installed upstream of the device. Further pre-treatment may be necessary depending on the feedwater quality.

If the total hardness is more than 10°dH and/or the ratio of temporary hardness to total hardness is more than 80%, a pre-treatment (e.g. softening system) is required to ensure the reliability and performance of the product.

4.1 Responsibility of the operator

- The installation and operating instructions must be kept in the immediate vicinity of the filter system and be accessible at all times.
- The filter system may only be operated in a technically perfect and operationally safe condition.
- The information in the installation and operating instructions must be followed in full.

4.2 Warranty and exclusion of liability

The information and recommendations given as well as the local drinking water regulations and disposal regulations applicable to the area of application must be observed. All information and notes in this installation and operating manual account for the applicable standards and regulations, state-of-the-art technology as well as our expertise in water treatment.

The filter cartridge comes with a 2-year warranty.

No liability is assumed for damages and consequential damages arising from:

- Failure to comply with the instructions in this installation and operating manual
- Improper use
- Improper, faulty installation
- Improper commissioning, operation or maintenance
- Mechanical damage to the device
- Unauthorised modifications
- Technical changes
- Use of non-approved components
- Neglecting to perform the required service and replacement tasks

4.3 Qualified personnel

Only trained persons and specialist personnel may install, commission and maintain the filter system.

- Trained personnel have been informed about the tasks assigned to them and the possible dangers of improper use and behaviour.
- Skilled personnel are able to install, commission and maintain the filter system as a result of their technical training, knowledge and experience as well as knowledge of the relevant regulations.

4.4 Pressure

For optimal functioning, the device requires a certain minimum operating pressure. In addition, the water pressure should not exceed the maximum permissible pressure.

⚠ CAUTION!

The feedwater pressure must always be between 0.15 and 0.4 MPa (1.5 and 4 bar) as measured directly at the RO unit.

If the pressure exceeds 0.4 MPa (4 bar), then a pressure reducing valve must be installed.

If the pressure drops below 0.15 MPa (1.5 bar), then a pressure booster system must be installed upstream.

- We recommend installing a stop valve at the infeed side of the device to enable the feedwater flow to be interrupted for service purposes.
- The on-site installation should use at least DN 10. Smaller feed pipes may result in operational stoppage due to insufficient water pressure or flow rates (e.g. when flushing the reverse osmosis membranes).
- Installing a pressure reducer can reduce the flow.

5 Installation and assembly

5.1 Requirements for installation

- Choose a place to set up the device that allows simple connection to the water supply network.
- There must also be a drain connection and a separate mains socket (220–240 V, 50–60 Hz) nearby.
- The device must be electrically connected to a grounded mains socket. The voltage supply and the requisite feedwater pressure must be guaranteed stable.
- Interference emissions (voltage peaks, high-frequency electromagnetic fields, spurious and voltage oscillations, etc.) of the surrounding electrical installation may not exceed the maximum values set in the standard EN 61000-6-4.

⚠ CAUTION!

Quality of the supply pipeline: All materials used in the permeate area must be corrosion-resistant.

Read *Technical data* [2] and *Operating and safety instructions* [4] prior to installation.

5.2 Unpacking the Planmeca ClinicPatrol

Remove the device from the packaging and check that the delivery is complete and undamaged.

⚠ CAUTION!

Defective parts must be replaced immediately.

Work hygienically.

5.3 Hydraulic installation

⚠ CAUTION!

Only hoses that comply with DVGW W 543 may be used to connect the device.

When installing accessories (hoses, connection sets), observe the installation dimensions and bending radii.

- The RO device is to be set up and operated with installed magnet plates.
- The hoses of the device are to be connected free of tension.
- The water connections must be connected in a watertight manner.
- The concentrate line is to be routed to the on-site wastewater connection with an air gap in accordance with EN 1717 and then connected there. The flexible hoses may not exhibit any constriction in cross section. During installation, ensure that the concentrate and permeate lines are correctly connected.
- The concentrate facing supplied must be connected to the concentrate outlet. The hose is then connected to the concentrate facing using the supplied 90° JG 8 mm fitting.

ⓘ NOTE!

Before using the RO unit, we recommend checking the water pre-treatment (e.g. in-house water softening system, central water processing of waterworks). This measure is necessary to improve the efficiency and service life of your RO membrane.

Please discard the permeate produced during the first 10 minutes after each new installation, initial commissioning or membrane change.

Reduction of the temperature by 1 °C results in a reduction of the permeate output of the membranes of approximately 3 percent.

Please also observe the provisions in the operating manual of the external pre-filter used.

As a general rule, we recommend using softened water to extend the service life and operational safety of the RO membrane.

6 Operating the reverse osmosis unit

6.1 Switching on the RO device

⚠ DANGER!

Always disconnect the voltage supply of the unit by removing the plug or fuse if hard wired during any maintenance and electrical work.

📌 NOTE!

The socket outlet shall be installed near the equipment and shall be easily accessible.

The device may only be operated with the supplied power cable.

- The RO device must be connected to the hoses and the electrical socket.
- Insert the mains plug (220–240 V/50–60 Hz).
- Open the cock for the feedwater supply.
- The RO device is switched on with the device switch on the back of the unit. A green LED indicates that the RO device is switched on (POWER ON).

6.2 Hygiene concept and stoppages

There are two hygiene concepts designed to ensure the maximum service life of the diaphragms:

- Stop delay: After each production run, the pump continues to run for approx. 10 s to rinse the diaphragm with tap water. This ensures that the conductivity on the concentrate side of the diaphragm falls back down to the input conductivity. This reduces stagnation peaks and prevents calcification of the membrane. The volume of waste water generated by this process is approx. 330 ml.
- Permeate rinsing: When the permeate conductivity exceeds a preset value, the permeate is returned via the membrane. If the value falls below the setpoint, demineralisation starts.
- Hygienic requirement: For use in clinical areas (e.g. in the dental area), an additional germ barrier must be installed after the device (e.g. an ultrafiltration device such as the BWT bestcare mini).

If the system is out of service without power for an extended period of time (for example, during holidays), it is necessary to flush the device for ten minutes. Ensure that the rinse water can drain freely. The flushing water must be discarded.

6.3 Removing/replacing the RO cartridge

6.3.1 Disassembling/assembling the Planmeca Patrol 5

- Switch off the device (on/off switch on the back of the device). Make sure that the status LED is no longer lit.
- Remove the new Planmeca Patrol 5 from the packaging and remove the hygiene cap.
- 6 ■ Before installing the Planmeca Patrol 5, write the date of installation and the replacement date (after 12 months at the latest) on the type plate of the filter cartridge.
- 8 ■ Tilt the device backwards slightly to have better access to the filter cartridge to be replaced.
- Turn the old filter cartridge clockwise to remove from the filter head.
- 9 ■ Screw the new filter cartridge anti-clockwise into the filter head.
- Power the unit back on and check the system for leaks.
- If the Planmeca Patrol 5 filter cartridge has been replaced, reset the post-filter counter in the app.

📌 NOTE!

10 After each Planmeca Patrol 5 replacement, the device must be rinsed for 5 minutes.

6.3.2 Disassembling/assembling the Planmeca Patrol 14 membrane

- Switch off the device (on/off switch on the back of the device). Make sure that the status LED is no longer lit.

- Remove the new Planmeca Patrol 14 membrane from the packaging and remove the hygiene cap.
- Before installing Planmeca Patrol 14 Membrane, write the date of installation and the replacement date (after 24 months at the latest) on the type plate of the filter cartridge.
- Tilt the device backwards slightly to have better access to the filter cartridge to be replaced.
- First, screw the Planmeca Patrol 5 clockwise out of the filter head.
- Screw the blind cap (bestmax blindcap 812162 / FS00Y30A00) into the free space for the Planmeca Patrol 5.
- You can now screw the Planmeca Patrol 14 membrane out of the device in a clockwise direction.
- Now screw the new Planmeca Patrol 14 membrane into the device in an anti-clockwise direction.
- The water meter (*membrane water counter*) for the Planmeca Patrol 14 membrane membrane cartridge can then be reset in the app. The Planmeca Patrol 14 membrane is now rinsed automatically.
- Once the rinsing process is complete, the blind cap has to be replaced again with the Planmeca Patrol 5.

6

7

8

9

10

11

ⓘ NOTE!

If the permeate volume flow decreases, the membrane element must be replaced. Furthermore, we recommend replacing the membrane every 12 months.

6.4 Installation and operation of the RO app

6.4.1 Installing the app

If you have not yet installed the BWT RO app on your mobile phone, scan the QR code below. This will take you to the website from which you can download the app. The app can also be downloaded from the following website: <https://www.bwt-group.net/bwt/download/bwt-RO-Systems.html>. **The app was developed only for iOS end devices.**

ⓘ NOTE!

Make sure that the Bluetooth connection of your end device is activated.

During the initial installation, after downloading the app, proceed as follows to activate the app: Settings / General / Device management / "BWT Aktiengesellschaft" / "Trust BWT Aktiengesellschaft"

The app is a purely offline application. No data is transmitted to BWT.

6.4.2 Using the app




The operating instructions of the BWT RO app can be found in the appendix.

7 Troubleshooting

7.1 Overview of the status and alert LED



Status	LED colour	Device state
Working	Lit up in green	Device in operation
Ready	Pulsates green	Device is in standby
Working	Lit up in yellow	WCF alert triggered
Working	Lit up in blue	Mobile device connected to RO via Bluetooth
Service	Lit up in red	Service required
Fault	Lit up in red or flashes red	Error message as displayed in the app (<i>Troubleshooting [7.2]</i>)

7.2 Troubleshooting

App screen	LED flashing mode	Cause	Measure
	 <p>Status LED lights up in red and a warning signal sounds.</p>	<p>Pressure drop in permeate outlet or pressure sensor faulty.</p>	<p>Check the permeate line for possible negative pressure / check the preload pressure in the pressure tank.</p> <p>Disconnect the device from the power and water supply.</p> <p>Contact your service technician.</p>
	 <p>Status LED lights up in red and a permanent warning signal sounds.</p>	<p>Water leakage inside the device or condensation</p>	<p>Disconnect the device from the power and water supply.</p> <p>Contact your service technician.</p>
	 <p>The status LED repeatedly flashes twice in succession</p>	<p>No flow or insufficient flow</p>	<p>Check the water pipe parameters.</p> <p>Check the inlet valve.</p> <p>Check that the pre-filter is not blocked.</p> <p>The device will restart automatically.</p>
	 <p>The status LED flashes.</p>	<p>Motor overheating protection has tripped.</p> <p>The pump may be broken.</p>	<p>The device starts by itself when the motor has cooled down.</p> <p>The device no longer starts by itself if the pump is broken.</p>

App screen	LED flashing mode	Cause	Measure
 <p>PLANMECA</p> <p>replace best demin</p> <p>0.1</p> <p>TRANSFER PUMP II</p> <p>ClinicPatrol</p>	<p>A warning tone sounds three times in succession. The status LED lights up in red.</p>	<p>Planmecca Patrol 5 has not been replaced; it is likely to need replacing. The conductivity no longer falls below the limit value.</p>	<p>Replace the Planmecca Patrol 5 and press the red symbol.</p>
 <p>PLANMECA</p> <p>transfer pump coding</p> <p>0.1</p> <p>TRANSFER PUMP II</p> <p>ClinicPatrol</p>	 <p>The status LED flashes 7 times in a row.</p>	<p>The transfer pump does not work or is not connected.</p>	<p>The device will start by itself when the transfer pump motor has cooled down.</p>
 <p>PLANMECA</p> <p>fault inlet h2o meter</p> <p>0.1</p> <p>TRANSFER PUMP II</p> <p>ClinicPatrol</p>	 <p>The status LED repeatedly flashes three times in succession.</p>	<p>The input water meter is blocked or has no power supply.</p>	<p>Permeate is still produced, but the volume of the input water is not recorded.</p> <p>Check the power supply.</p> <p>Contact your service technician.</p>
 <p>PLANMECA</p> <p>level sensor fault</p> <p>0.1</p> <p>TRANSFER PUMP II</p> <p>ClinicPatrol</p>	 <p>Status LED lights up in red and a warning signal sounds.</p>	<p>Problems with the tank fill level sensor.</p>	<p>Contact your service technician.</p>

2.7 ppm Patrol 5	14.2 °C	2.6 ppm MEMBRANE
0.0 bar OUTLET PRESS	108 DAYS	
fault liter MEMBRANE	0035 hours PUMP ON	
03657 liter Patrol 5	000 liter/hour OUTLET FLOW	
FAULT INLET H2O METER		

App screen	LED flashing mode	Cause	Measure
	 <p>The status LED flashes 9 times in a row.</p>	<p>The tank is overfilled. The fill level sensor indicates a value that is too high.</p>	<p>When the tap is opened, the vessel is emptied.</p>
	 <p>The status LED flashes 8 times in a row.</p>	<p>Problems with the membrane conductivity, it does not fall below the set threshold.</p>	<p>The device must be switched off to replace the membrane cartridge. Click the red area in the app to re-set.</p>
	 <p>The status LED flashes 10 times in a row.</p>	<p>Conductivity limit of the demineralisation filter exceeded. Planmecca Patrol 5 demineralisation filter cartridge is likely to need changing soon.</p>	<p>Replace the Planmecca Patrol 5 and press the red symbol.</p>

8 Repairs and maintenance

You have purchased a product that is durable and easy to service. To guarantee faultless functioning and optimal water quality, the unit must be maintained at regular intervals (at least once per year) by an authorised service technician.

If the product malfunctions during the warranty period, contact your contract partner or the installation company and quote the device type and serial number (see technical specifications or the type plate on the device).

ⓘ NOTE!

Before performing work on electrical components or opening the housing, it is mandatory that the mains plug be pulled and both the water infeed and the permeate line be closed off in order to guarantee that the unit is voltage-free.

Whenever maintenance is performed, the connection lines and the device must be checked for damage.

8.1 Wearing parts

Wearing parts must be replaced within the prescribed maintenance intervals.

Maintenance work	Responsibility	Interval
General visual inspection	Customer	Weekly
Leak-tightness inspection	Customer	Weekly
Cleaning with a damp cloth	Customer	As necessary
Conductivity (with external measuring device)	Customer/service	At least once per year
Replacement of the external pre-filter insert (particle filter [optionally available])	Customer/service	Depending on the pre-filter used
Replacement of the reverse osmosis filter cartridge	Customer/service	Once a year (recommended), after 2 years at the latest
Replacement of the demineralisation filter cartridge	Customer/service	As required, at the latest every 12 months
Replacement of filter head	Service	After 5 years, at the latest after 10 years

8.2 Cleaning

Clean your device with a damp cloth and a mild cleaning agent. To protect the surfaces of the device, do not use alcohols, bleach or solvents.

8.3 Standard IEC 60335-1

- This device is not intended to be operated by persons not in possession of full physical, sensory and mental faculties (including children) nor by persons without relevant experience or knowledge. Prior to using the device, personnel must be instructed in its use and given clear operating instructions by an expert specialist.
- The device is to be secured against access by children.
- To avoid hazards, in the event of any damage to the power cable, it must be replaced by the manufacturer, a service partner of the manufacturer or a similarly qualified person.
- Please visually inspect the water hoses for damage.

① NOTE!

In accordance with BVG A3 (VBG4), the unit must be tested for electrical safety every 4 years.

The reverse osmosis filter cartridge is subject to the "Pressure Equipment Directive" 2014/68/EU from 27 June 2014. The RO unit fulfils the requirements set forth in article 3, section 3, and was designed and manufactured according to the applicable engineering good practices.

The device does not have a CE mark in accordance with article 6, section 5, of the directive 2014/68/EU, but the attached CE Declaration of Conformity is valid.

8.4 Disposal

The device consists of various materials which need to be disposed of properly.

Please contact your contract partner for an expert and environmentally friendly disposal. Please do not dispose of depleted batteries in general household waste.

Any electronic parts should be disposed of only at authorised recycling centres (2012/19/EU). Observe the applicable national regulations on disposal of electric devices.

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EN: Select the device from the list. If multiple devices are present, note the ID number corresponding to the correct device.

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
EN: Once the Planmeca Clinic Patrol is selected, the screen below will be displayed. If the device is connected to the app, the status LED on the front of the device is BLUE.

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
EN: Once a connection to the device has been established, the screen below will be displayed. The device is ready for permeate production. The app shows "ready".

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
EN: During use, the app shows the status "working" with a blue spinning circle. Reverse osmosis is in operation.

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
EN: If you press the “working” button during operation, the following operating data appears: Conductivity of Planmeca Patrol 5 and Planmeca Patrol 14 membrane: If the target value has not yet been reached, the display flashes “ $\mu\text{S}/\text{cm}$ MEMBRANE” and no permeate is produced; Outlet pressure according to Planmeca Patrol 5; running time since commissioning; Planmeca Patrol 14 membrane and Planmeca Patrol 5 water meter; operating hours pump; water meter output

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EN: After each production run, the membrane is rinsed with inlet water for 10 seconds.

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
EN: A limit value can be set. When this limit value is reached, the device must be serviced. If the alarm is no longer applicable, it can be deactivated.

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
EN: Enter the following password to access the service menu: 05310.

If a value has been changed, the "save" button must be pressed, and an audible signal sounds. Press the "next" button to proceed to the next level. Press the "exit" button to leave the service menu and return to the operating data display.

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